

# Food & Animal Solutions



**SECURE** contactless.

# Food & Animal Solutions



## Safeguarding food, livestock and companion animals

As commerce becomes increasingly global, the task of tracking and managing animals and food products becomes ever-more critical and challenging. Animals cannot be stacked and stored like other products. They may be transported over long distances, often mingling with others from different locations. And farm animals, zoo animals and household pets alike are all potential carriers of diseases that can threaten animal and human health and lives.



*Secure contactless technology from HID can help identify lab animals, prevent the spread of disease and even track a pet's veterinary history.*

As with animals, tracking food products through processing and transport contributes to consumer health by identifying both product origins and possible problems in processing, handling, or transport. Aside from health concerns, tracking can also be critical for animals and food products that are particularly valuable or that require special protection. These can range from rare wines to family pets to breeding stock in zoos or on farms.

HID Global is a world leader in the development of radio frequency secure contactless technology used to manage and safeguard critical assets like these. HID Global provides a full range of specialized tags for use in managing farm, laboratory, exotic, and companion animals as well as high-value food products. Our innovative products meet all applicable standards and regulations and are fully interoperable with other

standardized components and systems.

### Meeting market demands

HID Global has long recognized consumers' concerns regarding food safety and quality. Farmers and processors have implemented rigorous systems for meeting consumer demands, but one of the biggest challenges has been tracking individual animals. This is critical for pinpointing animal origins, controlling the spread of disease, preventing fraud, and managing costs. In addition, national and regional bodies have implemented strict regulations for the handling of animals and animal-based food products.

Animals pass through many hands as they move from farm to exporter to auction house to processor. Because their movement cannot be precisely controlled, older visual technologies like tattoos and non-electronic ear-tags made tracking difficult. Secure contactless identification technologies supplement or augment these existing technologies, allowing animals to be tagged with radio frequency transponders and be identified even where direct line-of-sight is not practical.

### Establishing global standards

The effectiveness of contactless identification technology depends on establishment of, and compliance with, global standards. The memory coding structure and



communication protocols for contactless identification technology in the animal market are defined by ISO standards I 1784 and I 1785. As a member of the ISO committee responsible for these standards, HID Global has been in the forefront of their development and implementation and is committed to ongoing standardization of contactless technologies.

### **Non-food applications**

While food animals are the most obvious application for this technology, other applications are growing in importance. These include tagging zoo and laboratory animals as well as family pets. Secure identification of animals helps ensure the validity of research. It helps identify lost or stolen animals and aids in their return. It helps track animals' veterinary history, ensure proper care, and prevent the spread of disease. In the European Union, chips implanted under the skin of pets, coupled with an EU pet passport, permits travel of animals among EU member states without quarantine.

### **A world of options**

Effective contactless identification requires the right tools for each application. HID Global offers transponders designed to meet virtually any need including livestock animals of various sizes, pet identification, laboratory and exotic animals, and even fish. Transponders are incorporated into ear tags for combined visual/radio

frequency applications. Glass tags, sized from 4 x 34 mm for large animals like cattle to 2 x 10 mm for small animals like mice and snakes, can be safely and permanently injected subcutaneously or embedded in a ceramic bolus permanently placed into an animal's stomach. These ISO-compliant devices operate at 134.2 KHz and are available in read-only or one-time-programmable (OTP) configurations. Where appropriate, implanted transponders are coated to prevent migration.

HID also produces transponders for tracking food products—animal carcasses, cheeses, and perishable produce. In addition, HID makes contactless tags to track high-value bottled beverages and specialty items through the distribution chain.

### **Industry leadership**

HID Global is committed to the advancement of contactless automatic identification and data capture (AIDC) technologies. Our automated production lines and high-volume production capabilities are unmatched in the industry. With multiple ISO 9001:2000 certified and CC EAL 5+ certified manufacturing plants, we produce over 250 million units a year. Our R&D capabilities are second to none, and we work with end-users, developers and equipment suppliers, and system integrators to identify emerging needs and develop innovative solutions.



*Integrated contactless tracking systems can increase productivity through monitoring and nutritional controls, protect livestock from the threat of disease and ensure safety through precisely controlled processing and distribution.*

# Food & Animal Solutions

## Proven performance

HID Global has over 45 years of manufacturing experience and over 20 years of expertise in contactless technology development and manufacturing. Users and developers worldwide rely on HID Global for responsive solutions, fast turnaround, dependable quality, and exceptional performance.

## Partnerships and consulting services

At HID Global, we work in partnership with customers and system integrators to develop state-of-the-art systems to meet the most stringent real-world demands. Our components are fully compatible with the widest range of readers and back-end data systems and can, if necessary, be customized to meet special requirements.

## Animal Identification Selection Guide



	Livestock Identification				Pet and Exotic Animal Identification	
<b>Composition</b>	Glass tag	Glass tag	E-unit rod	E-unit disc	E-unit rod	Glass tag
<b>Tagging Method</b>	subcutaneous	ceramic bolus	ear tag	ear tag	customized implant solutions	subcutaneous
<b>Size</b>	2.12 mm x 12 mm 3.15 mm x 15.5 mm 4 mm x 22 mm 4 mm x 34 mm	4 mm x 34 mm	1.8 mm x 15 mm	27 mm 24 mm 28 mm & custom	1 mm x 8 mm	2.12 mm x 10 mm 2.12 mm x 12 mm

© 2009 HID Global. All rights reserved. HID and the HID logo are trademarks or registered trademarks of HID Global in the U.S. and/or other countries. All other trademarks, service marks, and product or service names are trademarks or registered trademarks of their respective owners.

mkt\_20apr2009\_f&a\_broch\_en



**SECURE contactless.**

[hidglobal.com](http://hidglobal.com)

**Corporate North America**  
15370 Barranca Pkwy  
Irvine, CA 92618  
U.S.A.  
Tel: 1 800 237 7769  
Tel: +1 949 732 2000  
Fax: +1 949 732 2360

**Contact for Food & Animal ID**  
HID Global Switzerland SA  
Rte Pra Charbon 27  
1614 Switzerland  
Tel: +41 908 011 00  
Fax: +41 21 908 01 03

**Asia Pacific**  
19/F 625 King's Road  
North Point  
Island East  
Hong Kong  
Tel: +852 3160 9800  
Fax: +852 3160 4809

**Latin America**  
Circunvalación Ote. #201 B  
Despacho 2  
Col. Jardines del Moral  
Leon 37160, Gto.  
Mexico  
Tel: +52 477 779 1492  
Fax: +52 477 779 1493

**Europe, Middle East & Africa**  
Haverhill Business Park  
Phoenix Road  
Haverhill, Suffolk  
CB9 7AE  
England  
Tel: +44 1440 714 850  
Fax: +44 1440 714 840